**QuikRead go<sup>®</sup> CRP Control Low**  
**QuikRead go<sup>®</sup> CRP Control**  
**QuikRead go<sup>®</sup> CRP Control High**

**Intended use**

QuikRead go CRP Control Low, QuikRead go CRP Control and QuikRead go CRP Control High are intended for quality control of CRP assays by the QuikRead go Instrument. **For *in vitro* diagnostic use.**

#### Principles of the procedure

Use the control in the same way as an unknown **plasma/serum** specimen would be used, following the instructions for the assay procedure.

#### Reagents

The control is human-based material with added color. It should be handled in the same way as patient sample. The control is stable for 15 minutes in the cuvette.

#### Preparation

The control is ready-to-use. Allow the material to stabilise for 30 minutes at room temperature (18...25°C). **Mix well.** Do not shake.

#### Storage

Stable until the expiry date when stored at 2...8°C. Opened vial: Stable for 2 months at 2...8°C. Do not use any vial that has visible evidence of microbial growth.

#### Expected values

The individual result should be within the range indicated on the control vial label. Each laboratory should establish its own mean and precision data.

|   |                          |                         |                       |
|---|--------------------------|-------------------------|-----------------------|
| <b>QuikRead go<sup>®</sup> CRP Control Low</b>  | <b>Kontrolle Niedrig</b> | <b>Kat. -Nr. 153765</b> | <b><i>Deutsch</i></b> |
| <b>QuikRead go<sup>®</sup> CRP Control</b>      | <b>Kontrolle</b>         | <b>Kat. -Nr. 153764</b> |                       |
| <b>QuikRead go<sup>®</sup> CRP Control High</b> | <b>Kontrolle Hoch</b>    | <b>Kat. -Nr. 153763</b> |                       |

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#### Warnings and precautions

- Reagents contain <0.1% sodium azide, which is not considered a harmful concentration. Sodium azide liberates toxic gas when in contact with acids. Azides may react with metal plumbing, forming explosive compounds. Azide build-up can be avoided by flushing with a large volume of water when disposing of the reagents.
- Avoid using the last drops of the control because opening of the vial many times may result in evaporation, which may bias the results. If the control is not in the acceptable range, take a new vial.
- Do not use the product after the expiry date marked on the outer package.
- Do not exceed the stability periods for opened reagents.
- The control contains material of human origin, which has been tested and found negative for human immunodeficiency virus antibodies (HIV 1 and 2), hepatitis C virus (HCV) and hepatitis B surface antigen (HBsAg). As a precaution, handle the control as if capable of transmitting infectious agents.

#### Disposal

- Disposal of contents according to national and local law.
- Materials of the components
  - Paper: Instructions for use
  - Cardboard: Outer box
  - Glass: Control vials
  - Plastic: Vial caps
- When used in accordance with Good Laboratory Practice, the good occupational hygiene and the instructions for use, the reagents supplied should not present a hazard to health.

|   |                          |                         |                       |
|---|--------------------------|-------------------------|-----------------------|
| <b>QuikRead go<sup>®</sup> CRP Control Low</b>  | <b>Kontrolle Niedrig</b> | <b>Kat. -Nr. 153765</b> | <b><i>Deutsch</i></b> |
| <b>QuikRead go<sup>®</sup> CRP Control</b>      | <b>Kontrolle</b>         | <b>Kat. -Nr. 153764</b> |                       |
| <b>QuikRead go<sup>®</sup> CRP Control High</b> | <b>Kontrolle Hoch</b>    | <b>Kat. -Nr. 153763</b> |                       |

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