

Line assays	Product number	Tests / Kit
RecomLine Borrelia IgM	130131 / 146954	20 and 200 tests
RecomLine Borrelia IgG	130132 / 146955	20 and 200 tests
CarL, fully automated strip processor	31094	1 pc

ELISA assays	Product number	Tests / Kit
Reagent C6 LYME IgM/IgG EIA	155181	96 tests
RecomWell Borrelia IgG	130186	96 tests
RecomWell Borrelia IgM	130187	96 tests
alphaWell Coxiella burnetii phase 1 IgG	151866	96 tests
alphaWell Coxiella burnetii phase 2 IgG	151867	96 tests
alphaWell Coxiella burnetii phase 2 IgM	151868	96 tests

Bead (Luminex) assays	Product number	Tests / Kit
recomBead Borrelia IgG 2.0	147419	96 tests
recomBead Borrelia IgM 2.0	147420	96 tests
recomBead CXCL13	151946	96 tests
Magpix instrument and Xponent 4.1 software	149969	1 pc

PCR assays	Product number	Tests / Kit
alphaCube TBE PCR	151961 / 151962	32 and 96 tests
alphaCube Borrelia	151969 / 151970	32 and 96 tests
alphaCube Babesia	151967 / 151968	32 and 96 tests
alphaCube Ehrlichia	151975 / 151976	32 and 96 tests
alphaCube Q-fever	151997 / 151998	32 and 96 tests

Rapid assays	Product number	Tests / Kit
ReaScan+ Rapid test Reader	154788	1 pc
ReaScan TBE IgM	147641	10 tests
ReaScan CXCL13	152375	10 tests
ReaScan+ C6 LYME IgG	154809	10 tests
ReaScan+ LYME IgM/IgG Combo	155180	10 tests

Your partner in detection of tick-borne infections



AIDIAN

Aidian Oy (Finland): info@aidian.fi / 050 380 9684 / For orders: tilaukset@aidian.fi
 Aidian Sweden AB: info@aidian.se / 08 623 04 00 / For orders: order@aidian.se
 Aidian Norway AS: info@aidian.no / 066 78 56 30 / For orders: ordre@aidian.no
 Aidian Denmark ApS: info@aidian.dk / 08 610 3020 / For orders: order@aidian.dk

9086-04EN Nordic, 06/2022

AIDIAN

Various tests to fit your needs in tick-borne disease diagnostics.

Ticks contain diverse pathogens and can pass these to humans. Our [wide range of products for tick-borne infections](#) support you with the correct diagnosis.

We offer several detection methods from high-throughput Borrelia screening to confirmation assays and easy to use rapid tests. Check out our portfolio below.

	Line	ELISA	Bead	PCR	Rapid test
Borrelia	x	x	x	x	x
CXCL13	-	-	x	-	x
TBE	-	-	-	x	x
Babesia	-	-	-	x	-
Coxiella burnetii	-	x	-	x	-
Ehrilica	-	-	-	x	-

Line

Strip-Immunoassay with specific recombinant antigens for the detection of IgG or IgM antibodies against Borrelia burgdorferi in human serum, plasma or CSF. The recomLine Borrelia uses immunodominant antigens leading to verification of early and late stages of Lyme disease. Fully automated CarL instrument available for easy sample-in-result-out processing of line tests.

ELISA

Enzyme immunoassays for detection of IgG or IgM antibodies from serum or CSF samples. ELISA tests can be used with open automated systems to allow minimum hands on time and fast results.

Bead

The advantages of screening and confirmation tests are combined in these Luminex® technology based borrelia IgG and IgM tests – specific antibody detection and quantification! Bead test also available for quantitative detection of CXCL13 in human cerebrospinal fluid.

PCR

Mikrogen alphaCube qualitative PCR assays for flexible analysis of various sample types. The PCR assays can be run on several open PCR systems and different extraction methods are suitable.

Rapid test

Reagent ReaScan®+ rapid tests for fast and easy antigen detection from serum or CSF with a small instrument that can be connected to your laboratory information system. ReaScan+ reader reports the result in numerical form allowing no room for subjective interpretations. Fast and simple!

